

# High-throughput picodroplet-based analysis of biosynthetic libraries

Carlos Bricio<sup>1</sup>, Xin Liu<sup>2</sup>, Xin Li<sup>2</sup>, Marian Rehak<sup>2</sup>, Tom Ellis<sup>3</sup> and Clive Smith<sup>2</sup>.








<sup>(1)</sup>Bioengineering, Imperial College London, London, United Kingdom, <sup>(2)</sup>Sphere Fluidics Limited, Cambridge, United Kingdom, <sup>(3)</sup>Centre for Synthetic Biology and Innovation, Imperial College London, London, United Kingdom

## Abstract

Biosynthesis of high value chemicals by engineered microorganisms is an application of synthetic biology that offers both economic and environmental advantages. This application is increasing the need for high-throughput screening tools that can facilitate the detection of the best performance among a library of designed microbes. For this reason we are developing a high-throughput, miniaturised Mass Spectrometry (MS) tool for profiling synthetic designed libraries. Combining microfluidics based picodroplet technology for cell encapsulation and sorting together with Mass Spectrometry we aim to rapidly screen, identify and retrieve the best cell "hits" among synthetic metabolic pathway libraries. Based on this novel approach we will be able to determine which construct has the genetic combination that gives the best biosynthesis performance. To test this new tool we have designed three libraries of two synthetic metabolic pathways using molecular engineering techniques<sup>(1,2,3)</sup>.

We chose two previously described synthetic pathways to produce non-natural amino acids<sup>(4, 5)</sup> and focus on improving their level of expression. Various strategies have been explored such as the use of homologous genes from other organisms, varying the DNA copy number, transcription levels or translation activity. Then, using a pioneering picodroplet-based technology<sup>(6)</sup> that enables not only the testing of up to 200,000 samples per day by MS, using miniaturised input volumes (400-700 pL), but also for retrieving identified "hits" in a reproducible manner, we will select single cells, analyse their production of this non-natural amino acids and finally select and recover the best performing clones among the different profiles obtained for further studies. This will enable new scientific breakthroughs, higher throughputs, lower screening costs, shorten design-build-test cycle and thus, be of interest to the current MS user base in the synthetic biology market and other sectors.

## DESIGN OF SYNTHETIC LIBRARIES



## Picodroplet MS Biochip

Nebulising N<sub>2</sub> gas

## Metal ESI-emitter




### MS screening

- MS Library Size**

We are using a PerkinElmer Axion 2 ESI mass spectrometer. One of the reasons for choosing this instrument is not only its high scan rate (100 Hz), but also its relatively high mass resolution (12000 FWHM at m/z 922), in order to be able to visualize individual picodroplets in the TIC. Ideally would like to have at least 10 scans to characterize each picodroplet, this leads to a MS screening throughput of up to 10 Hz (36,000 picodroplets/h).

### Modelling ESI MS Throughput

MS Screening Frequency (Hz)	10	5
Screening Time (Hours)	6	6
Library Size	216000	108000
Picodroplet #	5363457	2681729
Picodroplet generation (Hz)	1600	1600
Picodroplet Generation Time (Hours)	0.99	0.50
Picodroplet Diameter (μm)	98.5	98.5



## Conclusions & Future perspective

- It has been previously shown that microfluidic picodroplets can be used as reactors to study single bacterial proliferation<sup>(7,8)</sup>. The use of microfluidics in bacterial growth allow us to reduce the compartmentalization volume of bacterial cultures down to 500 pL. We have observed that picodroplet-based bacterial cultures do grow to a higher density than in 3 mL shake flasks. This technique will help to save screening costs and increase the number of samples that can be screened in a short period of time.
- We are currently increasing the variability of the synthetic libraries, optimising the microfluidics, growth media and developing the optics, barcoding strategy and software to analyse individual bacterial phenotypes that are producing non-natural amino acid.
- Ultimately, the recovery and analysis of the best hits after screening will allow to understand what combination of elements used in the synthetic library design is best for production and will help in further rounds of optimization.

## References

- Kittleson, J. T., Cheung, S., & Anderson, J. C. (2011). Rapid optimization of gene dosage in E. coli using DIAL strains. *J Biol Eng*, 5(1), 10.
- Engler, C., Gruetzner, R., Kandzia, R., & Marillonnet, S. (2009). Golden gate shuffling: a one-pot DNA shuffling method based on type IIs restriction enzymes. *PLoS One*, 4(5), e5553.
- Gibson, D. G., Young, L., Chuang, R.-Y., Venter, J. C., Hutchison, C. A., & Smith, H. O. (2009). Enzymatic assembly of DNA molecules up to several hundred kilobases. *Nature Methods*, 6(5), 343-9.
- Mehl, R. A., Anderson, J. C., Santoro, S. W., Wang, L., Martin, A. B., King, D. S., et al. (2003). Generation of a bacterium with a 21 amino acid genetic code. *J. Am. Chem. Soc.*, 125(4), 935-9.
- Zhang, K., Li, H., Cho, K. M., & Liao, J. C. (2010). Expanding metabolism for total biosynthesis of the non-natural amino acid L-homoalanine. *Proceedings of the Proc Natl Acad Sci USA*, 107(14), 6234-9.
- CA Smith, X Li, TH Mize, TD Sharpe, El Graziani, C Abell, WT Huck (2013). Sensitive, High Throughput Detection of Proteins in Individual, Surfactant-Stabilized Picoliter Droplets Using Nano-electrospray Ionization Mass Spectrometry. *Anal. Chem.* (85), 3812-3816.
- Martin, K., Henkel, T., Baier, V., Grodrian, A., Schön, T., Roth, M., et al. (2003). Generation of larger numbers of separated microbial populations by cultivation in segmented-flow microdevices. *Lab Chip* 3: 202-7
- Huebner, A., Srisa-Att, M., Holt, D., Abell, C., Hoffelder, F., deMello, A.J., and Edel, J.B. (2007) Quantitative detection of protein expression in single cells using droplet microfluidics. *Chem Commun (Camb)* 1218-20.